

Leedo[®]**HS-2205
Rotary Microtome****Features:**

- 1. It have clipper -built case, the user can put blade, wax block, or other consumable on the top of case. it is easy.
- 2. It do not need to lubricate and maintain the microtome frequently.
- 3. The hand wheel in the right can be locked in any position, and the handle can be locked in highest position, it is double safe for changing blade or tissue block.
- 4. The specimen Orienting system: axes X and Y 8°, axis 180°.
- 5. It adopt the special blade holder, it is more safe and convenient to use the disposable blade directly. It also can to longer the user life of disposable blade.

Main Technical Data:

- 1. Slice thickness range: 0μm ~ 60μm
- 2. Section thickness setting: 0-2 μm in 0.5 μm -increments
 - a)2-10μm in 1 μm-increments
 - b)10-20 μm in 2 μm-increments
 - c)20-60μm in 5 μm-increments
- 3. Horizontal specimen stroke: 20 mm
- 4. Vertical specimen stroke: 50mm
- 5. Slip precision: ± 5%
- 6. Maximal slice section: 50 × 45mm
- 7. Specimen adjustment direction: horizontal 8°, vertical 8°, rotating 360°
- 8. Movement range of blade holder base: front and back ±25mm, left and right ±25mm
- 9. Blade holder angle 0~10°

Net Weight: 30kg

Dimensions: 500*480*280mm

Accessories:**Standard Setting List**

- | | |
|----------------------------|-------|
| 1. Disposable Blade Holder | 1 pc |
| 2. Specimen Cassette Clamp | 1 pc |
| 3. Hexagonal Wrench | 1 set |
| 4. Brush | 1 pc |
| 5. Optional Accessories | |
| 6. Paraffin Wax Holder | |
| 7. Feather Microtome Blade | |

Package Info:

Each instrument packed Wooden Case
 Dimensions: 520 *520*450mm,
 Gross Weight: 41kg