



FOR DENTIST'S EXCELLENCE



KLT-6210 | N2+ Dental Unit

Hygiene Multifunctional
Excellent performance





LED Sensor Operation Light

Touchless on/off, stepless intensity control, prevent cross infection
Three axis rotation to meet different operation request
Detachable handle, easy for clean and sterilize
Light intensity: 7000-32000Lux
Stable illuminant, no RF interference, eye protection



Dentist Table

Wide side (650*300mm) table, meet the needs of dentist
Cleanable plastic table pad protect the table, longer life
Integral table handle more convenient and high efficient

Swiveling Instruments holder

Handpiece holder with 5 positions, pre-position for scaler or electric motor



Rotatable Water Unit Box

Provide more space for assistant working closely with dentist and patient at all time
Built-in tissue box, convenient for patients
Rotatable tray for storing patients personal effects
With two water bottle, sterilizing water supply system



Spittoon

Rotable & detachable ceramic spittoon (180 degree rotatable & Easy for clean)
Programable cup filler and spittoon rinsing system
Constant temperature warm water system
Pure water supply system



Detachable Suction Filter

Prevent bacteria from accumulating water inside with filter net, easy to take out for cleaning (With nylon filter, easy clean, corrosion resistance)



Patient Seat

One starter for air/water/electricity, ease your daily job
Polyether water & air tube, high quality, 98A strength
24V noiseless chair motor, with aluminum shell, longer life



PU·upholstery^{Shade guide}



Green-0#



Light Blue-1#



Middle Blue-3#



Dark Blue-4#



Mango Yellow-6#



Orange-7#



Claret-8#



Apple Green-9#



Black-10#



Red-11#



Pink-12#



Silver-13#



Grey-14#



Rose-15#



Nacarar-16#



Jewel Green-18#

FOSHAN ROSON MEDICAL INSTRUMENTS CO.,LTD

CONTACT US

Add: NO.9 Henggui Mid-Road, Lianhe, Industrial Zone, Luocun, Nanhai, Foshan, Guangdong Province.

Tel: +86-757-86407037

Fax: +86-757-86407038

Web: www.fsrosen.com

E-mail: sales001@fsrosen.com

