

YA-D5-1

Automatic Hospital ICU Bed



*Better care
Better comfort!*

YA-D5-1

Automatic Hospital ICU Bed

DESCRIPTION

MEDIK YA-D5-1 electric hospital bed is full electric operation type ICU bed to take care of the patient in ICU area. It will provide convenient support for both caregiver and patient. The road to recovery is often a critical and sensitive experience. Both patient and medical staff will have to work hand-in-hand to achieve their goal of nursing the patient back to health.

FEATURES

- Bed Frame Electrophoresis and powder coating
- Bed platform removable ABS platform
- Motor Quiet and robust electric actuators provide reliable operation
- Build-in controller can control all electrically operated functions and movements
- 4 ABS side rails with Tuck-away model, collapsible below the bed surface
- Head and Footboard ABS material, detachable
- Central brake system with a wheel diameter of 125mm
- Hand remote control
- Patient control inner constructed in side rails
- Nurse control located in the footboard, with side rail control and a functions lockout feature
- Auto-counter to provide the patient from sliding onto foot end
- Electric CPR function control by one single button on the side rails
- Manual CPR

STANDARD SPECIFICATION

Overall size	2200*980 *450-750mm
Lying size	1950*850 mm
Backrest adjustable	0-70°
Footrest adjustable	0-40°
Height adjustable	450-750mm
Trendelenburg	0-12°
Reverse-Trendelenburg	0-12°
Loading capacity	220kg



BED BOARD

The bed surface is detachable and has a handle for easy cleaning also it with vent holes for ventilation.



SIDE RAIL LOCK

Tuck-away full length railing meets IEC60601-2-52 safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.



SIDERAIL CONTROL

All functions can be operated by the built-in side rail control buttons, which are located inside and outside of the side rail, which is convenient for patients and medical staff to operate.



Foot board control panel

Nurse control panel embedded onto the foot board, railing control operation can be locked accordingly when position operation from patient is not wanted.



Detachable head and foot board

The detachable head and foot of the bed are equipped with safety locks to lock the head and foot board when it is not in use, and it can be quickly disassembled, which is convenient for medical staff to operate.



CENTRAL LOCK

Central lock system is convenient for the nurse to push the patient to the required position and quickly lock it with one key, which can appropriately reduce the workload of the medical staff.

